

RETAIL PRICE LIST



SiMPore Inc.  
150 Lucius Gordon Drive  
Suite 119  
West Henrietta, NY 14586

Tel: 1-888-323-NANO  
Tel: 1-585-214-0585  
Fax: 1-888-249-2935  
Email: sales@temwindows.com

**Domestic Shipping**

FedEx 2-Day \$25  
FedEx Overnight \$40

**International Shipping**

FedEx International Priority  
Please call for quote (typically \$65-\$100)

<b>Subtotal</b>	
<b>Sales Tax</b>	
<b>Shipping</b>	
<b>Total</b>	

**Ship To Address:**

Contact Person
Organization
Street
City,State Zip
Country
Phone #

**Bill To Address:**

Contact Person
Organization
Street
City,State Zip
Country

**Payment Information:**

**(Please do not fax or email credit card info, Call us with this information)**

VISA                       MASTERCARD

Name on Card:

CC Number:

Expiration:

3-Digit Code:

Purchase Order Number:

Payment Terms: Net 30, 1% late fee each month

RETAIL PRICE LIST

**PURE SILICON (amorphous or nanoporous)**

Part Number	Frame Thickness	Film Thickness	Description	RETAIL Price per 10	Packages of 10
US100-A05Q00	100 µm	5 nm	One Window: (1) 25x25 µm	<b>\$290</b>	
US100-A05Q33A	100 µm	5 nm	Nine Windows (SMALL): (8) 50x50, (1) 50x100 µm	<b>\$290</b>	
US100-A05L	100 µm	5 nm	Two Windows: (2) 50x1500 µm	<b>\$290</b>	
US100-A09Q33	100 µm	9 nm	Nine Windows: (8) 100x100 µm, (1) 100x350 µm	<b>\$250</b>	
US100-A09L	100 µm	9 nm	Two Windows: (2) 100x1500 µm	<b>\$250</b>	
US100-A15Q33	100 µm	15 nm	Nine Windows: (8) 100x100 µm, (1) 100x350 µm	<b>\$210</b>	
US100-A15L	100 µm	15 nm	Two Windows: (2) 100x1500 µm	<b>\$210</b>	
US100-P30Q05	100 µm	30 nm	<b>NANOPOROUS</b> , One Window: (1) 500x500 µm	<b>\$250</b>	
US100-P30Q33	100 µm	30 nm	<b>NANOPOROUS</b> , Nine Windows: (8) 100x100, (1) 100x350 µm	<b>\$250</b>	

RETAIL PRICE LIST

**SILICON NITRIDE (amorphous)**

Part Number	Frame Thickness	Film Thickness	Description	RETAIL Price per 10	Packages of 10
SN100-MP2Q05	100 µm	50 nm	<b>MICROPOROUS (2.0 µm pores on 1:1 pitch)</b> One Window: (1) 500x500 µm	<b>\$290</b>	
SN100-A05Q00	100 µm	5 nm	One Window: (1) 25x25 µm	<b>\$290</b>	
SN100-A05Q33A	100 µm	5 nm	Nine Windows (SMALL): (8) 50x50, (1) 50x100 µm	<b>\$290</b>	
SN100-A05L	100 µm	5 nm	Two Windows: (2) 50x1500 µm	<b>\$290</b>	PENDING
SN100-A10Q33	100 µm	10 nm	Nine Windows: (8) 100x100 µm, (1) 100x350 µm	<b>\$150</b>	
SN100-A20Q33	100 µm	20 nm	Nine Windows: (8) 100x100 µm, (1) 100x350 µm	<b>\$120</b>	
SN100-A20Q05	100 µm	20 nm	One Window: (1) 500x500 µm	<b>\$120</b>	
SN100-A50Q33	100 µm	50 nm	Nine Windows: (9) 100x100 µm	<b>\$90</b>	
SN100-A50Q10	100 µm	50 nm	One Window: (1) 1000x1000 µm	<b>\$90</b>	
SN100-A50Q05	100 µm	50 nm	One Window: (1) 500x500 µm	<b>\$90</b>	
SN100-A50Q01	100 µm	50 nm	One Window: (1) 100x100 µm	<b>\$90</b>	

RETAIL PRICE LIST

**SILICON DIOXIDE (amorphous)**

<b>Part Number</b>	<b>Frame Thickness</b>	<b>Film Thickness</b>	<b>Description</b>	<b>RETAIL Price per 10</b>	<b>Packages of 10</b>
SO100-A20Q33	100 $\mu\text{m}$	20 nm	Nine Windows: (9) 100x100 $\mu\text{m}$	<b>\$200</b>	
SO100-A40Q33	100 $\mu\text{m}$	40 nm	Nine Windows: (9) 100x100 $\mu\text{m}$	<b>\$150</b>	

## RETAIL PRICE LIST

### TEMVu™ Slides

SiMPore is pleased to offer standard and custom TEMVu™ slides made to order. Our retail price is \$20 per slide in addition to the base price of the TEM windows.

**For custom combinations - please contact us. Lead times are typically 3 business days.**

<b>TEMVu™ Part Number</b>	<b>RETAIL price per 10</b>	<b>Packages of 10</b>
TEMVu-SN100-A20Q05	<b>\$320</b>	
TEMVu-SN100-MP2Q05	<b>\$490</b>	

### ACCESSORIES

<b>Part Number</b>	<b>Description</b>	<b>RETAIL Price</b>	<b>Quantity</b>
GELBOX	Gelbox with 5x5 Grid	<b>\$5</b>	
K6TWZR	Plastic K6 TEM Window Handling Tweezer	<b>\$14</b>	